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ENGLISH

INSTRUCTION FOR USE

Name of MD: Wound drainage system (REDON)

BASIC UDI-DI: 859419979REDONXH

Model:

- Vacuum bottle for wound drainage Redon with tube (OR System)
- Vacuum bottle for wound drainage Redon without tube (Ward System)

PRODUCT DESCRIPTION:

The Wound drainage system (REDON) is a high-vacuum suction drainage system. It is a disposable medical device. The identical product is supplied in two versions: Sterile version and Non-sterile version. The product is made of polymer material. The product is pre-evacuated with an absolute maximum vacuum of 40 mbar (relative vacuum max. 900 mbar = 90 kPa). The system comes in 3 variations of drainage bottles 400 ml round, 200 ml round and 200 ml square, and a connecting tube with connectors.

- **Vacuum bottle for Redon wound drainage with tube (OR system)** consists of a 400ml or 200ml drainage bottle and a connecting tube (Large Lock connecting tube). The connecting tube is Ch18 with a Large-lock connector and Ch6-18 drain connector.
- **Vacuum bottle for Redon wound drainage without tube (Ward system)** (replaceable bottle) consists of a 400ml or 200ml drainage bottle. The connecting tube with a Large Lock connector is not part of this product.

The connecting tube (Large Lock) is a flexible transparent Ch18 tube with an outer diameter of 6 mm and a length of 1.2 m. It is made of soft natural PVC. The tube is equipped with a Large-Lock positive connector and a drain connector for connecting it to the drainage bottle to the suction drain and the tube clip. The connecting tube is connected to the drainage bottle with a Large-Lock connector. The connecting tube is connected to the drains through the drain connector.

The drain connector is a universal connector for connecting drains of various diameters (from Ch6 size to Ch18 size). The drain connector is indicated with the direction of suction and the Ch size of the individual drains. The terminations are glued to the connecting tube. The tube clip is designed for the possibility to close the flow in the connecting tube when applying it or when changing the drainage bottle.

Large Lock connecting tube:

Tube termination:	Designed to be connected with:	Second tube termination:	Designed to be connected with:
Drain connector CH 6 – CH18	Suction drain CH 6 – CH18 (It is not part of the product)	Large- Lock positive connector	Drainage bottle with Large- Lock negative connector

Size:

- 400 ml – round drainage bottle
- 200 ml – round drainage bottle
- 200 ml – square drainage bottle

The system is a single-use device. The device is suitable for use in MRI.

The sterile version of the products and its blister packaging are sterilised with ethylene oxide (EO).

Products supplied non-sterile are intended for use in sets and kits, and are supplied without blister packaging. Sterilisation with EO is performed after they are packed into a set. In this case, sterilisation is not included in our delivery.

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Intended use:

A source of vacuum energy associated with the removal of the extract into a drainage bottle.

Clinical benefit: The devices do not have a direct clinical benefit. An indirect clinical benefit is in connection with the drain. Devices help the drains with which they are used achieve their intended use.

Indications:

Not specified

Contraindications:

Not specified

Undesirable side effects:

Not associated with the use of the product

Intended patient population:

Global population of all age groups, including young people.

Intended users:

Surgeons and nurses that are qualified and familiar with the surgical techniques.

Warnings:

- The medical device is for single use. It is intended for one patient only. Reuse of a single-use device creates a potential risk of infection for the patient or user.
- Do not use if the packaging has been opened or damaged or the product is damaged.
- Do not use it after the end of its shelf life.
- Sterile unless the packaging has been opened or damaged.
- Check the blue vacuum indicator on the drainage bottle before use. If it is straight, the drainage bottle does not have a vacuum. Discard the drainage bottle and replace it with a spare drainage bottle.
- The medical device is not intended for the removal of body fluids or storage for subsequent return to the human body

Storage

- Store in the original packaging in a clean, dry, ventilated area away from direct light.
- To be stored in its packaging at temperatures between 5°C and 30°C.
- The storage temperature should not exceed 40°C for a long time (max. 30 days).
- During transport, the temperature should not be less than -20°C and should not exceed +50°C.

Directions for use

1. Check the integrity and completeness of the packaging.
2. Remove the drainage bottle from its packaging shortly before use.
3. Check the integrity and completeness of the product.
4. Use sterile equipment to take the drainage bottle out.
5. Insert the drain.
6. Cut the drain close to the needle at a 45° angle for easy connection the drain connector.
7. Cut the drain connector at the appropriate drain size.
8. Introduce the wound drain into the single drain connector, check the correct fixation.
9. Connect the connecting tube with the Large-Lock in the drainage bottle (Caution: loose connections can lead to a loss of vacuum). For connecting use proportionate force using fingers. Do not use instruments as they could damage the parts.
10. Open the sliding clamp on the connection piece of the drainage bottle. Check the suction function. If the vacuum is not spontaneously released from the drainage bottle after opening the white sliding clamp, squeeze the blue rubber part transversely. If the vacuum cannot be released, replace the drainage bottle.
11. Ensure an airtight wound seal.
12. Ensure that the vacuum indicator is completely compressed, which indicates that maximum negative pressure is available.

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13. Check the suction function.

14. The system is now put into operation and can be attached to the bed, in any position, by using the bed strap.

Connection of replacement drainage bottle


1. Replacement is necessary, e.g. when the vacuum indicator is completely expanded, which indicates minimum negative pressure available. Close the sliding clamp on the connecting tube.
2. Seal off the drainage bottle by closing the sliding clamp on the blue connection piece.
3. Disconnect the connecting tube by unscrewing the Large Lock.
4. Discard the full system, properly sealed, according to normal hospital policy as hazardous waste.
5. Take a replacement drainage bottle and repeat the steps from direction for use (points 9–14).
6. Open the sliding clamp on the connecting tube.

Disconnection

1. After the surgeon (or any other authorised person) has taken the decision to stop the wound drainage, close the tube clip on the connecting tube and the closure clip on the drainage bottle.
2. Before removing the drain from the patient's side, ensure that the pressure difference inside and outside the patient will be eliminated by letting the wound rest for 10–20 minutes.
3. Repeat steps 3 and 4 of the connection of a replacement drainage bottle.

Waste disposal: After disconnecting, close the drainage bottle and connecting tube and place them in a hazardous bio-waste collection container. Dispose of them as hazardous bio-waste in accordance with waste management regulations.

Symbols used on labelling:

Do not use if the packaging is damaged

Note: Any serious incident that has occurred in relation to the medical device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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